

**Test Format:** Triple Antibody Sandwich (TAS) ELISA

## Intended Use:

---

This capture antibody is intended for use with its corresponding detection antibodies (ECA).

Test performance using other detection antibodies has not been evaluated and is not guaranteed.

## Includes:

- Capture Antibody
- 96-well microtiter plates, solid

## India Contact:

### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222 Mobile: +91-9810521400 , Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Web: [www.lifetechindia.com](http://www.lifetechindia.com)