

Enzyme Conjugate for Johnsongrass mosaic virus (JGMV)



Test Label: Alkaline Phosphatase
Test Format: Triple Antibody Sandwich (TAS) ELISA

Intended Use:

This detection antibody (bottle A) and enzyme conjugate (bottle B) are intended for use with their corresponding capture antibody (CAB).

Test performance using other capture antiserum has not been evaluated and is not guaranteed.

Includes:

- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- [User Guide](#)

Option

Enzyme conjugate for Johnsongrass mosaic virus (JGMV), alkphos, 500 testwells
ECA 18300/0500

Enzyme conjugate for Johnsongrass mosaic virus (JGMV), alkphos, 1000 testwells
ECA 18300/1000

Enzyme conjugate for Johnsongrass mosaic virus (JGMV), alkphos, 5000 testwells
ECA 18300/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com