

Enzyme Conjugate for Banana bunchy top virus (BBTV)



[Click to view a larger image](#)

Test Label: Alkaline Phosphatase
Test Format: Triple Antibody Sandwich (TAS) ELISA

Intended Use:

This detection antibody (bottle A) and enzyme conjugate (bottle B) are intended for use with their corresponding capture antibody (CAB).

Test performance using other capture antiserum has not been evaluated and is not guaranteed.

Includes:

- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- [User Guide](#)

Option

Enzyme conjugate for Banana bunchy top virus (BBTV), alkphos, 500 testwells
ECA 24700/0500

Enzyme conjugate for Banana bunchy top virus (BBTV), alkphos, 1000 testwells
ECA 24700/1000

Enzyme conjugate for Banana bunchy top virus (BBTV), alkphos, 5000 testwells
ECA 24700/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com