

## Enzyme Conjugate for Tobacco etch virus (TEV)



**Test Label:** Alkaline Phosphatase

**Test Format:** Double Antibody Sandwich (DAS) ELISA

### Intended Use:

This enzyme conjugate is intended for use with its corresponding capture antibody (CAB).

Test performance using other capture antiserum has not been evaluated and is not guaranteed.

### Includes:

- Alkaline phosphatase enzyme conjugate
- [User Guide](#)

### Option

Enzyme conjugate for Tobacco etch virus (TEV), alkphos, 500 testwells  
ECA 49501/0500

Enzyme conjugate for Tobacco etch virus (TEV), alkphos, 1000 testwells  
ECA 49501/1000

Enzyme conjugate for Tobacco etch virus (TEV), alkphos, 5000 testwells  
ECA 49501/5000

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)