

## Enzyme Conjugate for *Pantoea stewartii* (Pstew)



**Test Label:** Alkaline Phosphatase  
**Test Format:** Triple Antibody Sandwich (TAS) ELISA

### Intended Use:

This detection antibody (bottle A) and enzyme conjugate (bottle B) are intended for use with their corresponding capture antibody (CAB).

Test performance using other capture antiserum has not been evaluated and is not guaranteed.

### Includes:

- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- [User Guide](#)

<input type="checkbox"/>	Option
<input type="checkbox"/>	Enzyme conjugate for <i>Pantoea stewartii</i> (Pstew), alkphos, 500 testwells ECA 52000/0500
<input type="checkbox"/>	Enzyme conjugate for <i>Pantoea stewartii</i> (Pstew), alkphos, 1000 testwells ECA 52000/1000
<input type="checkbox"/>	Enzyme conjugate for <i>Pantoea stewartii</i> (Pstew), alkphos, 5000 testwells ECA 52000/5000

### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)