

## Enzyme Conjugate for Iris yellow spot virus (IYSV)



**Test Label:** Alkaline Phosphatase  
**Test Format:** Double Antibody Sandwich (DAS) ELISA

### Intended Use:

This enzyme conjugate is intended for use with its corresponding capture antibody (CAB).

Test performance using other capture antiserum has not been evaluated and is not guaranteed.

### Includes:

- Alkaline phosphatase enzyme conjugate
- [User Guide](#)

### Option

Enzyme conjugate for Iris yellow spot virus (IYSV), alkphos, 500 testwells  
ECA 60500/0500

Enzyme conjugate for Iris yellow spot virus (IYSV), alkphos, 1000 testwells  
ECA 60500/1000

Enzyme conjugate for Iris yellow spot virus (IYSV), alkphos, 5000 testwells  
ECA 60500/5000

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)