

Anti-animal Enzyme Conjugate for Xanthomonas albilineans (Xalb)



Test Label: Alkaline Phosphatase
Test Format: Antigen coated plate (ACP) ELISA

Intended Use:

This anti-animal enzyme conjugate is intended for use with its corresponding detection antibody in ACP ELISA format.

Test performance using other capture antiserum has not been evaluated and is not guaranteed.

Includes:

- Alkaline phosphatase enzyme conjugate
- [User Guide](#)

<input type="checkbox"/>	Option
<input type="checkbox"/>	Anti-animal enzyme conjugate for Xanthomonas albilineans (Xalb), alkphos, 500 testwells ECA 71100/0500
<input type="checkbox"/>	Anti-animal enzyme conjugate for Xanthomonas albilineans (Xalb), alkphos, 1000 testwells ECA 71100/1000
<input type="checkbox"/>	Anti-animal enzyme conjugate for Xanthomonas albilineans (Xalb), alkphos, 5000 testwells ECA 71100/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com