

Enzyme Conjugate for Cymbidium mosaic virus (CymMV)



Test Label: Alkaline Phosphatase
Test Format: Triple Antibody Sandwich (TAS) ELISA

Intended Use:

This detection antibody (bottle A) and enzyme conjugate (bottle B) are intended for use with their corresponding capture antibody (CAB).

Test performance using other capture antiserum has not been evaluated and is not guaranteed.

Includes:

- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- [User Guide](#)

Option

CymMV enzyme conjugate, alkphos, 500 testwells
ECA 89501/0500

CymMV enzyme conjugate, alkphos, 1000 testwells
ECA 89501/1000

CymMV enzyme conjugate, alkphos, 5000 testwells
ECA 89501/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com