


Lipoamido-PEG12-acid

Cat. No.:	HY-141334
CAS No.:	2407442-47-9
Molecular Formula:	C ₃₅ H ₆₇ NO ₁₅ S ₂
Molecular Weight:	806.03
Target:	PROTAC Linker
Pathway:	PROTAC
Storage:	Please store the product under the recommended conditions in the Certificate of Analysis.



BIOLOGICAL ACTIVITY

Description	Lipoamido-PEG12-acid is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs ^[1] .
IC₅₀ & Target	PEGs
In Vitro	PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins ^[1] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.

REFERENCES

[1]. An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. *EBioMedicine*. 2018 Oct;36:553-562

Caution: Product has not been fully validated for medical applications. For research use only.