

ImmunoStrip® for Xanthomonas (Xan)

The *Xanthomonas* (Xan) ImmunoStrip® is used to detect the presence of Xanthomonas to the genus level and cannot differentiate species. It can detect Xanthomonas in fruits, ornamentals, and vegetables. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

Option
ImmunoStrip® for Xan, 5 strips and buffer filled bags ISK 14600/0005
ImmunoStrip® for Xan, 25 strips and buffer filled bags ISK 14600/0025
ImmunoStrip® for Xan ***BUFFER FILLED BAGS NOT INCLUDED*** STX 14600/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitamputa, Delhi – 110034 (India). Ph: +91-11-42208000, 42208111, 42208222 Mobile: 91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com