

ImmunoStrip® for Impatiens necrotic spot virus (INSV)

The *Impatiens necrotic spot virus* (INSV) ImmunoStrip® is used to detect the presence of INSV in ornamental and vegetable crops. INSV is a member of the Orthospovirus genus known for their enveloped, spherical-shaped virus particles. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

<input type="checkbox"/>	Option
<input type="checkbox"/>	ImmunoStrip® for Impatiens necrotic spot virus (INSV), 5 strips & buffer filled bags ISK 20501/0005
<input type="checkbox"/>	ImmunoStrip® for Impatiens necrotic spot virus (INSV), 25 strips & buffer filled bags ISK 20501/0025
<input type="checkbox"/>	ImmunoStrips® only for Impatiens necrotic spot virus (INSV) ***BUFFER FILLED BAGS NOT INCLUDED*** STX 20501/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com