

## Agdia ImmunoStrip® for Tobacco streak virus (TSV)

The *Tobacco streak virus* (TSV) ImmunoStrip® is used to detect the presence of TSV in ornamental and vegetable crops. TSV is a member of the Ilarvirus genus known for their non-enveloped, quasi-spherical-shaped virus particles. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

### *Included in ISK Kit:*

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

### *Included in STX Item:*

- ImmunoStrips®
- [User Guide](#)

Option
ImmunoStrip® for Tobacco streak virus (TSV), 5 strips & buffer filled bags  ISK 25500/0005
ImmunoStrip® for Tobacco streak virus (TSV), 25 strips & buffer filled bags  ISK 25500/0025
ImmunoStrips® only for Tobacco streak virus (TSV) ***BUFFER FILLED BAGS NOT INCLUDED***  STX 25500/0025

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (India). Ph: +91-11-42208000, 42208111, 42208222 Mobile: 91-11-42208444  
 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)