

ImmunoStrip® for Squash mosaic virus (SqMV)



Test Label: NA
Test Format: ImmunoStrip®
Capture Reagent: Polyclonal
Detection Reagent: Polyclonal

Intended Use:

The SqMV ImmunoStrip® is an on-site tool to quickly identify SqMV-infected plants with or without external disease symptoms.

The test requires no special equipment or expertise, making it ideal for use in the field.

Kit (ISK) Includes:

- ImmunoStrips®
- Sample bags containing SEB1 buffer
- [User Guide](#)

ImmunoStrips® Only (STX) Includes:

- ImmunoStrips®
- SEB1 extraction buffer purchased separately

Option
ImmunoStrip® for Squash mosaic virus (SqMV), 5 strips & buffer filled bags ISK 26400/0005
ImmunoStrip® for Squash mosaic virus (SqMV), 25 strips & buffer filled bags ISK 26400/0025
ImmunoStrips® only for Squash mosaic virus (SqMV) ***BUFFER FILLED BAGS NOT INCLUDED*** STX 26400/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com