

ImmunoStrip® for Pepper mild mottle virus (PMMoV)

The *Pepper mild mottle virus* (PMMoV) ImmunoStrip® is used to detect the presence of PMMoV in pepper leaves and seeds. This ImmunoStrips® has been shown to reliably detect 1 PMMoV-infected seed in a subsample of 249 PMMoV-free pepper seeds. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

Option
ImmunoStrip® for Pepper mild mottle virus (PMMoV), 5 strips and buffer filled bags ISK 37502/0005
ImmunoStrip® for Pepper mild mottle virus (PMMoV), 25 strips and buffer filled bags ISK 37502/0025
ImmunoStrip® for Pepper mild mottle virus (PMMoV), 25 strips ONLY *** BUFFER FILLED BAGS NOT INCLUDED *** STX 37502/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (India). Ph: +91-11-42208000, 42208111, 42208222 Mobile: 91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com