

ImmunoComb® for CMV, INSV, TMV, TSWV



Test Label: NA
Test Format: ImmunoCombs®
Capture Reagent: NA
Detection Reagent: NA

Intended Use:

The ImmunoComb® joins together the most requested Agdia tests for plant virus detection.

The grower can test for *Cucumber mosaic virus* (CMV), *Tomato spotted wilt virus* (TSWV), *Impatiens necrotic spot virus* (INSV) and *Tobacco mosaic virus* (TMV) all at one time.

Samples are ground in a special bag that allows for easy extraction. The comb of four test strips is then easily placed into the sample. The test takes just minutes to form clear yes or no results.

Kit (ISK) Includes:

- 24 ImmunoCombs®
- Sample bags containing SEB1 buffer
- [User Guide](#)

ImmunoCombs® Only (STX) Includes:

- 24 ImmunoCombs®
- SEB1 extraction buffer purchased separately

Option

ImmunoComb® Kit for CMV, INSV, TMV, & TSWV
ISK 37800/0024

ImmunoCombs® only for CMV, INSV, TMV, TSWV ***BUFFER FILLED
BAGS NOT INCLUDED***
STX 37800/0024

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com