

Overview

Synonyms	Histone H3, human;
Description	Recombinant protein, full-length, expressed in <i>E. coli</i> and purified using FPLC. Purity 90% by SDS-PAGE and Commassie blue staining. Accession No. NP_003520.
Species	Human
Source	<i>E. coli</i>
Sequence	MARTKQTARK STGGKAPRKQ LATKAARKSA PATGGVKKPH RYRPGTVALR EIRRYQKSTE LLIRKLPFQR LVREIAQDFK TDLRFQSSAV MALQEACEAY LVGLFEDTNL CAIHAKRVTI MPKDIQLARR IRGERA

Properties

Measured Molecular Weight	Theoretical M.W.=15.4 kDa
Formulation	100 µg of human histone H3 quantitated using the molar absorption coefficient (276 nm) of 4640. Lyophilized powder; reconstitute in sterile water or buffer of choice.
Storage	Stable for 2 years at 4°C from date of shipment. Once rehydrated, aliquot and store at -20°C.

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com