

Overview

Synonyms	Hhg1
Description	<p>Sonic Hedgehog (Shh) is a member of the Hedgehog (Hh) family of highly conserved proteins which are widely represented throughout the animal kingdom. In mammal, there are three related Hh proteins, Sonic (Shh), Desert (Dhh) and Indian (Ihh). They share a high degree of amino-acid sequence identity (e.g., Shh and Ihh are 93% identical). Sonic Hedgehog plays a role in cell growth, cell specialization, and the normal shaping (patterning) of the body. Shh is also important for development of the brain and spinal cord (central nervous system), eyes, limbs, and many other parts of the body.</p> <p>Recombinant mouse Sonic Hedgehog (C25II) (rmShh) produced in <i>E. coli</i> is a single non-glycosylated polypeptide chain containing 175 amino acids. A fully biologically active molecule, rmShh has a molecular mass of 19.7 kDa analyzed by reducing SDS-PAGE and is obtained by proprietary chromatographic techniques.</p>
Accession No	Q62226
Species	Mouse
Source	<i>E. coli</i>
Biological Activity	ED ₅₀ < 2.0 µg/ml, measured by its ability to induce alkaline phosphatase production by C3H/10T1/2 (CCL-226) Cells, corresponding to a specific activity of > 500 units/mg.
Sequence	<p>IIGPGRGFGK RRHPKLTPL AYKQFIPNVA EKTLGASGRY EGKITRNSER FKELTPNYNP DIIFKDEENT GADRLMTQRC KDKLNALAIS VMNQWPGVKL RVTEGWDEDG HHSEESLHYE GRAVDITTSR RDRSKYGMLA RLAVEAGFDW VYYESKAHII CSVKAENSVV AKSGG</p>

Properties

Measured Molecular Weight	19.7 kDa, observed by reducing SDS-PAGE.
Purity	> 95% by SDS-PAGE and HPLC analyses.
Formulation	Lyophilized after extensive dialysis against PBS.
Reconstitution	Reconstituted in ddH ₂ O at 100 µg/ml.
Endotoxin Level	< 0.2 EU/µg, determined by LAL method.
Storage	Lyophilized recombinant mouse Sonic Hedgehog (C25II) (rmShh) remains stable up to 6 months at lower than -70°C from date of receipt. Upon reconstitution, rmShh should be stable up to 2 weeks at 4°C or up to 3 months at -20°C.
Note	For research use only

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com