

Overview

Synonyms	Lymphopoietin 1(LP-1), pre-B cell factor
Description	Interleukin-7 (IL-7) , also known as lymphopoietin 1 and pre-B cell factor, is a hematopoietic growth factor belonging to the IL-7/IL-9 family. It is produced by keratinocytes, dendritic cells, hepatocytes, neurons and epithelial cells. IL-7 binds and signals through IL-7 receptor, a heterodimer consisting of IL-7 receptor alpha and common gamma chain receptor. IL-7 plays a role in regulating early B cell and T cell development. It is also important for optimal dendritic cell-T cell interaction.
Accession No	P13232
Species	Human
Source	CHO
Biological Activity	ED ₅₀ <0.2ng /ml, measured in a cell proliferation assay using 2E8 cells.
Sequence	DCDIEGKDGK QYESVLMVSI DQLLDSMKEI GSNCLNNEFN FFKRHICDAN KEGMFLFRAA RKLRLQFLKMN STGDFDLHLL KVSEGTILL NCTGQVKGRK PAALGEAQPT KSLEENKSLK EQKKLNDLCF LKRLQLQEIKT CWNKILMGTK EHHHHHHH

Properties

Measured Molecular Weight	24-34 kDa, observed by reducing SDS-PAGE.
Purity	> 95% as analyzed by SDS-PAGE.
Formulation	Lyophilized after extensive dialysis against PBS.
Reconstitution	Reconstituted in ddH ₂ O or PBS at 100 µg/ml.
Endotoxin Level	< 0.2 EU/µg, determined by LAL method.
Storage	Lyophilized recombinant human Interleukin-7 (IL-7), His remains stable up to 6 months at lower than -70°C from date of receipt. Upon reconstitution, human Interleukin-7 (IL-7), His should be stable up to 1 week at 4°C or up to 2 months at -20°C.
Note	For research use only

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com