

Overview

Synonyms	SDF-1 ² , Stromal-Cell Derived Factor-1, CXCL12, PBSF
Description	<p>SDF-1[±] and SDF-1², members of the chemokine[±] subfamily that lack the ELR domain, were initially identified using the signal sequence trap cloning strategy from a mouse bone-marrow stromal cell line. SDF-1[±] and SDF-1² cDNAs encode precursor proteins of 89 and 93 amino acid residues, respectively. Both SDF-1[±] and SDF-1² are encoded by a single gene and arise by alternative splicing. The two proteins are identical except for the four amino acid residues that are present in the carboxy-terminus of SDF-1² and absent from SDF-1[±]. SDF-1/PBSF is highly conserved between species, with only one amino acid substitution between the mature human and mouse proteins. SDF-1/PBSF acts via the chemokine receptor CXCR4 and has been shown to be a chemoattractant for T-lymphocytes, monocytes, pro- and pre-B cells, but not neutrophils. Mice lacking SDF-1 or CXCR4 have been found to have impaired B-lymphopoiesis, myelopoiesis, vascular development, cardiogenesis and abnormal neuronal cell migration and patterning in the central nervous system.</p> <p>Recombinant Mouse SDF-1²/CXCL12 produced in <i>CHO</i> cells is a polypeptide chain containing 78 amino acids. A fully biologically active molecule, rm SDF-1²/CXCL12 has a molecular mass of 8.5 kDa analyzed by reducing SDS-PAGE and is obtained by chromatographic techniques.</p>
Source	<i>CHO</i>
Biological Activity	The EC ₅₀ value of mouse SDF-1β/CXCL12 on Ca ²⁺ mobilization assay in CHO-K1/G±15/mCXCR4 cells (human G±15 and mCXCR4 stably expressed in CHO-K1 cells) is less than 2.5 µg/ml.
Sequence	Lys ²² -Met ⁹³ (accession #: P40224)

Properties

Measured Molecular Weight	8.5 kDa, observed by reducing SDS-PAGE.
Purity	> 95% as analyzed by SDS-PAGE.
Formulation	Lyophilized after extensive dialysis against PBS.
Reconstitution	Reconstituted in ddH ₂ O or PBS at 100 µg/ml.
Endotoxin Level	< 0.2 EU/µg, determined by LAL method.
Storage	Lyophilized recombinant Mouse SDF-1²/CXCL12 remains stable up to 6 months at lower than -70°C from date of receipt. Upon reconstitution, Mouse SDF-1 ² /CXCL12 should be stable up to 1 week at 4°C or up to 3 months at -20°C.
Note	For research use only

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com