

Life Technologies (India) has in public interest decided to provide Oxygen Concentrators to those who need it at a no-profit basis in light of the Covid19 pandemic. The equipment shall be transferred at cost price on declaration by the purchaser that it shall not be resold/used in the market for any kind of profit activity. Kindly note that very limited machines are available with us of the model JAY-5AW. We would like to support the fellow citizens to the best of our abilities. Please contact us at customerservice@lifetechindia.com along with the Covid19 positive report or proof of your healthcare centre. If we could be of any assistance at this hour of crisis, we would be most obliged.

MODEL: JAY-3AW/JAY-5AW



Low noise,
takes up less space



USA imported humidifier,
6PSI pressure release



Innovative big
LCD display

Optional

1. Nebulizer(Atomization)
2. Low Purity alarm
3. Low purity alarm with low flow alarm
4. High temperature alarm
5. Pulse oximeter
6. Maintenance alarm
7. Reset Switch
8. Remote control monitoring system (WIFI/Bluetooth)

Specification

Flow Rate	0-3LPM / 0-5LPM
Purity	93% (±3%)
Outlet Pressure	0.04-0.07MPA(6-10PSI)
Sound Level	< 38db
Power Consumption	< 300w / < 350w
LCD Display	Switch times, operating pressure, present working time, accumulated working time, presetting time from 10mins to 40hours
Alarm	Power Failure alarm, High&low pressure alarm
Net Weight	14Kgs / 15Kgs
Size	348×280×510 mm

Machine manufactured by Longfian Scitech Co. Ltd.

Model: JAY-5AW

Voltage and Frequency: 230V/50HZ

Power cord: EU plug

Machine Label and user manual: Included

Certificate

The Certification Body of
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization
LONGFIAN SCITECH CO., LTD.
2F&3F, East Section, Building 12
Power valley pioneer park
No. 369 Huiyang street
Baoding
071051 Hebei
China

has established and applies a quality management system for medical devices
for the following scope:

Manufacture and Distribution of Medical Oxygen Concentrators

Proof has been furnished that the requirements specified in

EN ISO 13485:2016

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2019-04-18
Certificate Registration No.: SX 60136154 0001
An audit was performed. Report No.: 16802952 007
This Certificate is valid until: 2022-03-18

Certification Body



Date 2019-04-18



TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
Tel.: +49 221 806-1371 Fax: +49 221 806-3935 e-mail: cert-validity@de.tuv.com <http://www.tuv.com/safety>



TÜVRheinland®

EC Certificate
Directive 93/42/EEC Annex V
Production Quality Assurance
Medical Devices

Registration No.: DD 60136153 0001

Report No.: 16802952 009

Manufacturer:

LONGFIAN SCITECH CO., LTD.
2F&3F, East Section, Building 12
Power valley pioneer park
No. 369 Huiyang street
Baoding
071051 Hebei
China

Products:

Medical Oxygen Concentrators

Replaces Approval, Registration No.: DD 60108301 0001

Expiry Date:

2024-03-18

The Notified Body hereby declares that the requirements of Annex V of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex V, section 4 of the aforementioned directive. For placing on the market of class IIb and class III devices covered by this certificate an EC type-examination certificate according to Annex III is required.

Effective Date:

2019-03-19

Date:

2019-03-18



TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg

TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.