

## Product Description

Adipose-Derived Stem Cell Attachment Medium is designed to promote adipose-derived stem cells in stromal vascular fraction attachment during cell culture. This medium is fully supplemented with fetal bovine serum, growth factors. This medium is sterile filtered and is ready for use.

<b>Catalog No.:</b>	MDADAM
<b>Package Size:</b>	250 mL, 500 mL
<b>Bioassay:</b>	Each lot was tested for the ability to promote attachment of adipose-derived stem cells in culture
<b>Sterility test:</b>	Negative for bacteria, yeast, fungi and mycoplasma.
<b>Storage &amp; stability:</b>	Avoid exposure to light. Suitable for 6 months at 4°C.

India Contact:

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura,  
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)

Web: [www.lifetechindia.com](http://www.lifetechindia.com)