

Preadipocyte Growth Medium is designed to promote preadipocytes attachment, spreading and proliferation in culture. This medium is fully supplemented with fetal bovine serum, growth factors. This medium is sterile filtered and is ready for use.

**Catalog No.:** MD-0004-500ML

**Package Size:** 500 mL

**Bioassay:** Each lot was tested for the ability to promote plating and proliferation of preadipocytes in culture.

**Sterility test:** Negative for bacteria, yeast, fungi and mycoplasma.

**Storage & stability:** Avoid exposure to light. Suitable for 6 months at 4°C.

India Contact:

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura,  
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)

Web: [www.lifetechindia.com](http://www.lifetechindia.com)