

## Product Description

Osteogenic Differentiation Medium is designed to efficiently induce adipose-derived stem cells differentiation into osteoblasts. This medium is fully supplemented and containing 10% fetal bovine serum. This medium is sterile filtered and is ready for use.

- Catalog No.:** MD-0006-250ML  
**Package Size:** 250 mL  
**Bioassay:** Each lot was tested for the ability to induce differentiation of adipose-derived stem cells into osteoblasts with accumulation of extracellular calcium deposit within 2~3 weeks.  
**Sterility test:** Negative for bacteria, yeast, fungi and mycoplasma.  
**Storage & stability:** Avoid exposure to light. Suitable for 6 months at 4°C.

India Contact:

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura,  
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)

Web: [www.lifetechindia.com](http://www.lifetechindia.com)