

Human Neural Stem Cell Growth Medium

SKU: MD-0024

SKU	Package Size
MD-0024-100ML	100 mL
MD-0024-500ML	500 mL

Description

Human Neural Stem Cell Growth Medium is designed to efficiently maintain human neural stem cells (hNSCs) in both suspension culture and adherent culture. This medium is fully supplemented, sterile filtered and is ready for use.

Catalog No.:	MD-0024
Package Size:	100 mL, 500 mL
Bioassay:	Each lot was tested for the ability to maintain NSCs in adherent culture for 3 passages.
Product Use:	Human Neural Stem Cell Growth Medium is for research use only. It is not approved for human or animal use, or for application in in vitro diagnostic procedures.
Sterility test:	Negative for bacteria, yeast, fungi and mycoplasma
Storage & stability:	Avoid exposure to light. Store at -20°C for up to 12 months. Once thawed, the medium can be stored at 4°C for up to 2 weeks.

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com