

## Product Description

Pericyte Growth Medium is designed to promote pericyte proliferation in culture. This medium is fully supplemented with fetal bovine serum, growth factors and antibiotic. This medium is sterile filtered and is ready for use.

<b>Catalog No.:</b>	MD-0030-500ML
<b>Package Size:</b>	500 mL
<b>Bioassay:</b>	Each lot was tested for the ability to promote proliferation of pericyte in culture
<b>Sterility test:</b>	Negative for bacteria, yeast, fungi and mycoplasma
<b>Storage &amp; stability:</b>	When stored in the dark at 4°C, the medium is stable for 6 months. It is recommended to store small aliquots at -20°C for longer time (1-2 years).

India Contact:

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura,  
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)

Web: [www.lifetechindia.com](http://www.lifetechindia.com)