

Product Description

Trabecular Meshwork Cell Growth Medium is designed to promote trabecular meshwork cells proliferation *in vitro*. This medium is fully supplemented with fetal bovine serum, growth factors without antibiotic. This medium is sterile filtered and is ready for use.

Catalog No.:	MD-0057
Quantity:	500 mL
Bioassay:	Each lot was tested for the ability to promote proliferation of trabecular meshwork cells in vitro.
Product Use:	Trabecular Meshwork Cell Growth Medium is for research use only. It is not approved for human or animal use, or for application in in vitro diagnostic procedures.
Sterility test:	Negative for bacteria, yeast, fungi and mycoplasma.
Storage & stability:	Suitable for 1 year at 4°C.

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com