

## Product Description

Sertoli Cell Growth Medium is designed to promote Sertoli cell proliferation in culture. This medium is fully supplemented with growth factors and antibiotics. This medium is sterile filtered and is ready for use.

**Catalog No.:** MD-0091

**Quantity:** 500 mL

**Bioassay:** Each lot was tested for the ability to promote proliferation of Sertoli cells in culture.

**Sterility test:** Negative for bacteria, yeast, fungi and mycoplasma.

**Storage & stability:** When stored in the dark at 4°C, the medium is stable for 6 months. It is recommended to store small aliquots at -20°C for longer time (>1 years).

India Contact:

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura,  
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)

Web: [www.lifetechindia.com](http://www.lifetechindia.com)