

Product Description

Human Dopaminergic Neuron Maturation Medium is designed for *in vitro* culture of iPSC-derived dopaminergic neurons. This medium is fully supplemented. This medium is sterile filtered and is ready for use.

Catalog No.:	MD-0105
Package Size:	100 mL
Bioassay:	Each lot was tested for the ability to culture human iPSC-derived dopaminergic neurons.
Product Use:	Human Dopaminergic Neuron Maturation Medium is for research use only. It is not approved for human or animal use, or for application in <i>in vitro</i> diagnostic procedures.
Sterility test:	Negative for bacteria, yeast, fungi and mycoplasma
Storage & stability:	Avoid exposure to light. Stored at -20°C for one year. Stored for 2 weeks at 4°C.

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com