

CspB ELISA Kit

Test Label: Peroxidase

Test Format: Double Antibody Sandwich (DAS) ELISA

Intended Use:



This ELISA was developed for the qualitative detection of the CspB protein in corn leaf. This ELISA is validated to detect the CspB protein expressed from Monsanto's transgenic corn event MON87460.

When testing leaf, Agdia recommends testing coleoptile leaf tissue when possible. Lower OD values are expected when testing leaf material beyond the coleoptile stage.

This test shows no cross-reaction with Bt-Cry1Ab, Bt-Cry1A.105, Bt-Cry2A, Bt-Cry3Bb1, NPTII, or CP4 EPSPS transgenic proteins.

Includes:

- Antibody-coated 96-well microtiter plates
- Peroxidase enzyme conjugate, concentrated
- RUB6 conjugate diluent
- TMB substrate solution
- Positive control
- PBST buffer, 20X concentrate
- [User Guide](#)

Option
CspB ELISA Kit, 480 testwells PSP 35500/0480
CspB ELISA Kit, 4800 testwells PSP 35500/4800

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com