

eCry3.1Ab ELISA Kit



Test Label: Peroxidase

Test Format: Double Antibody Sandwich (DAS) ELISA

Intended Use:

This ELISA was developed to detect the presence or the absence of the eCry3.1Ab protein expressed in corn seed, leaf or grain. This ELISA is validated to detect the eCry3.1Ab protein expressed in products containing Syngenta's corn event 5307 (Agrisure Duracade®).

This test shows no cross-reaction with the following genetically modified traits/events in available corn seed GA21 (mEPSPS), Bt11 (Cry1Ab), mCry3A,

or MIR162 (Vip3A).

Includes:

- Antibody-coated 96-well microtiter plates
- Peroxidase enzyme conjugate, concentrated
- RUB6 conjugate diluent
- TMB substrate solution
- Positive control
- PBST buffer, powder
- **User Guide**

Option

eCry3.1Ab Elisa Kit, 480 testwells
PSP 49000/0480

eCry3.1Ab bulk Elisa Kit, 4800 testwells
PSP 49000/4800

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com