

## Bt-Cry1Ac & Bt-Cry2A ELISA Kit

**Part Number:** PSM 14900/4800

- Overview

**Test Label:** Alkaline Phosphatase & Peroxidase

**Test Format:** Double Antibody Sandwich (DAS)  
ELISA

**Capture Reagent:** NA

**Detection Reagent:** NA



### Intended Use:

This dual-trait ELISA is validated to detect the presence or absence of Bt-Cry1Ac & Bt-Cry2A proteins expressed in transgenic cotton seed and leaf. The test does not distinguish between Bt-Cry1Ab and Bt-Cry1Ac proteins.

This test shows no cross-reaction with Bt-Cry1F, CP4 EPSPS, PAT/bar, Vip3A, Bt-Cry1C or NPTII transgenic proteins.

### Includes:

- Antibody-coated 96-well microtiter plates
- Alkaline phosphatase and peroxidase enzyme conjugate mix, concentrated
- RUB6 enzyme conjugate diluent
- pNPP substrate solution, 1X ready to use
- TMB substrate solution
- Positive control
- PBST buffer, powder
- [User Guide](#)

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)