

Bt-Cry34Ab1 ELISA Kit

Test Label: Perox

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: NA

Detection Reagent: NA

Intended Use:

This qualitative Bt-Cry34Ab1 ELISA is designed for single seed and single leaf detection of the Cry34Ab1 protein expressed in Herculex[®] RW and Herculex[®] XTRA seed corn. This double antibody sandwich ELISA was developed and validated by Agdia and approved for testing Herculex[®] traits by DowAgro Sciences.

Agdia GMO Trait ELISA's are commonly used by plant breeders, universities and private testing labs throughout the world.

Herculex[®], Herculex[®] RW and Herculex[®] XTRA are registered trademarks of Dow AgroSciences, LLC

Includes:

- Antibody-coated 96-well microtiter plates
- Peroxidase enzyme conjugate, concentrated
- RUB6 enzyme conjugate diluent
- TMB substrate solution
- Positive control
- PBST buffer, 20X concentrate
- [User Guide](#)



Bt-Cry34Ab1 ELISA Kit, 480 testwells,PSP 04500/0480

Bt-Cry34Ab1 ELISA Kit, 4800 testwells,PSP 04500/4800

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com