

Bt-Cry2A ELISA Kit

Test Label: Perox

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: NA

Detection Reagent: NA

Intended Use:

This ELISA was developed for the qualitative detection of the Bt-Cry2A protein in corn leaf or seed. This ELISA is validated to detect the Bt-Cry2Ab protein expressed from Monsanto's transgenic corn event, MON89034.

The Bt-Cry2A antibody recognizes Bt-Cry2A protein and shows no cross-reaction with Bt-Cry3Bb1, Bt-Cry34Ab, Bt-Cry1Ab, Bt-Cry1Ac, Bt-Cry1F, CP4 EPSPS (Roundup Ready[®], GA21, NPT II, mBt-Cry3A (Agrisure[®] RW) or PAT (LibertyLink[®]).

LibertyLink[®] is a registered trademark of Bayer CropScience

Roundup Ready[®] is a registered trademark of Monsanto

Agrisure[®] RW is a registered trademark of Syngenta Group

Includes:

- Antibody-coated 96-well microtiter plates
- Peroxidase enzyme conjugate, concentrated
- RUB6 conjugate diluent
- TMB substrate solution
- Positive control
- PBST buffer, 20X concentrate
- [User Guide](#)



India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com

Bt-Cry2A ELISA Kit, 96 testwells,PSP 05801/0096

Bt-Cry2A ELISA Kit, 288 testwells,PSP 05801/0288

Bt-Cry2A ELISA Kit, 480 testwells,PSP 05801/0480

Bt-Cry2A ELISA Kit, 4800 testwells,PSP 05801/4800

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com