

mBt-Cry3A ELISA Kit

Part Number:PSP 06700/4800

• [Overview](#)

TestLabel:Peroxidase

Test Format:Double Antibody Sandwich (DAS)
ELISA

Capture Reagent:N/A

Detection Reagent:N/A



Intended Use:

The Agdia mBt-Cry3A ELISA is a qualitative immunoassay which effectively detects the modified Cry3A protein (mCry3A) expressed from Agrisure[®] RW (MIR604) seed corn and leaves. The mCry3A protein can also be detected from the hybrid corn stacked-trait varieties, Agrisure[®] CB/LL/RW and Agrisure[®] GT/RW.

The mBt-Cry3a antibody recognizes mBt-Cry3A protein and shows no cross-reaction with Bt-Cry3Bb1, Bt-Cry34Ab, Bt-Cry1Ab, Bt-Cry1Ac, Bt-Cry1F, Bt-Cry9C, CP4 EPSPS (Roundup Ready[®]) or PAT (LibertyLink[®]).

Agdia GMO Trait ELISA's are complete kits and are commonly used by plant breeders, universities and private testing labs throughout the world.

LibertyLink[®] is a registered trademark of Bayer CropScience

Roundup Ready[®] is a registered trademark of Monsanto Company

Agrisure[®] RW is a registered trademark of Syngenta Group Company

Includes:

- Antibody-coated 96-well microtiter plates
- Peroxidase enzyme conjugate, concentrated
- RUB6 enzyme conjugate diluent
- TMB substrate solution
- Positive control
- PBST buffer, powder
- [User Guide](#)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com