

## RoundUp Ready<sup>®</sup> ELISA Kit

**Test Label:** Perox

**Test Format:** Double Antibody Sandwich (DAS) ELISA

**Capture Reagent:** NA

**Detection Reagent:** NA

### Intended Use:

The Agdia Roundup Ready<sup>®</sup> ELISA is designed for determining the presence or absence of the CP4 EPSPS protein expressed in Roundup Ready<sup>®</sup> corn, cotton and soybeans. This ELISA is validated for single seed and single leaf detection or bulk composite testing (0.1% for all crops).

Please see product insert for additional information.

### Includes:

- Antibody-coated 96-well microtiter plates
- Peroxidase enzyme conjugate
- TMB substrate solution
- Positive control
- Nonfat dried milk
- PBST buffer, 20X concentrate
- Tween 20
- [User Guide](#)



RoundUp Ready<sup>®</sup> ELISA Kit, 96 testwells,PSP 74000/0096

RoundUp Ready<sup>®</sup> ELISA Kit, 288 testwells,PSP 74000/0288

RoundUp Ready<sup>®</sup> ELISA Kit, 480 testwells,PSP 74000/0480

RoundUp Ready<sup>®</sup> ELISA Kit, 4800 testwells,PSP 74000/4800

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)