

Vip3A ELISA Kit

Test Label: Perox

Test Format: Double Antibody Sandwich (DAS) ELISA

Intended Use:

This ELISA was developed for the qualitative detection of the Vip3A protein in corn leaf or seed. This ELISA is validated to detect the Vip3A protein expressed in commercial products containing Syngenta's transgenic corn event MIR162 (Agrisure Viptera[®]).

This test shows no cross-reaction with Bt-Cry1Ab, Bt-Cry1A.105, Bt-Cry2A, mBt-Cry3A, Bt-Cry1F, Bt-Cry3Bb1, Bt-Cry34/35Ab1, ecry3.1Ab, CspB, amy797E, PAT/pat, NPTII, GA21 or CP4 EPSPS transgenic proteins.

Includes:

- Antibody-coated 96-well microtiter plates
- Peroxidase enzyme conjugate, concentrated
- RUB6 conjugate diluent
- TMB substrate solution
- Positive control
- PBST buffer, 20X concentrate
- [User Guide](#)



Vip3A ELISA Kit, 480 testwells,PSP 83500/0480

Vip3A ELISA Kit, 4800 testwells,PSP 83500/4800

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com