

Pan Path

REMBRANDT®

Locus specific identifiers ISH detection assay

LSIs are important for identifying specific genes or genetic markers that may be involved in diseases. ISH analysis using LSIs is a powerful tool for studying genetic rearrangements and chromosomal abnormalities. The REMBRANDT® Locus specific identifiers detection assays are able to detect and localize genes.

The detection of aneuploidy with locus specific identifier ISH probes has important clinical applications in various fields and can provide valuable diagnostic and prognostic information.

Visit www.panpath.nl for more information or contact your local distributor

Product list*

LOCUS SPECIFIC IDENTIFIERS;COMPLETE KITS

CAT NR.

- REMBRANDT® LSI 3q26.2-FISH orange detection assay C730K.3000.YY**
- REMBRANDT® LSI 3q26.2-ISH detection assay C730K.XXXX.YY**

- REMBRANDT® LSI 3q26.33-FISH orange detection assay C731K.3000.YY**
- REMBRANDT® LSI 3q26.33-ISH detection assay C731K.XXXX.YY**

- REMBRANDT® LSI 3pter-FISH orange detection assay C732K.3000.YY**
- REMBRANDT® LSI 3pter-ISH detection assay C732K.XXXX.YY**

- REMBRANDT® LSI 3qter-FISH orange detection assay C733K.3000.YY**
- REMBRANDT® LSI 3qter-ISH detection assay C733K.XXXX.YY**

- REMBRANDT® LSI 4p16.3-FISH orange detection assay C734K.3000.YY**
- REMBRANDT® LSI 4p16.3-ISH detection assay C734K.XXXX.YY**

*All products are also available as stand-alone RUO products

**Available in 5 or 10 tests

LOCUS SPECIFIC IDENTIFIERS;COMPLETE KITS**CAT NR.**

- REMBRANDT® LSI 7p11.2-FISH orange detection assay **Bestseller!** C736K.3000.YY**
- REMBRANDT® LSI 7p11.2-ISH detection assay C736K.XXXX.YY**

- REMBRANDT® LSI 9p21-FISH orange detection assay C739K.3000.YY**
- REMBRANDT® LSI 9p21-ISH detection assay C739K.XXXX.YY**

- REMBRANDT® LSI 10p14-FISH orange detection assay C740K.3000.YY**
- REMBRANDT® LSI 10p14-ISH detection assay C740K.XXXX.YY**

- REMBRANDT® LSI 12q15-FISH orange detection assay **Bestseller!** C741K.3000.YY**
- REMBRANDT® LSI 12q15-ISH detection assay C741K.XXXX.YY**

- REMBRANDT® LSI 13q12-FISH orange detection assay C742K.3000.YY**
- REMBRANDT® LSI 13q12-ISH detection assay C742K.XXXX.YY**

- REMBRANDT® LSI 15q11.2-FISH orange detection assay C743K.3000.YY**
- REMBRANDT® LSI 15q11.2-ISH detection assay C743K.XXXX.YY**

*All products are also available as stand-alone RUO products

**Available in 5 or 10 tests

LOCUS SPECIFIC IDENTIFIERS;COMPLETE KITS**CAT NR.**

- REMBRANDT® LSI 17q12-FISH orange detection assay **Bestseller!** C747K.3000.YY**
- REMBRANDT® LSI 17q12-ISH detection assay C747K.XXXX.YY**

- REMBRANDT® LSI 22q11.21-FISH orange detection assay C748K.3000.YY**
- REMBRANDT® LSI 22q11.21-ISH detection assay C748K.XXXX.YY**

- REMBRANDT® LSI Xp22.31-FISH orange detection assay C750K.3000.YY**
- REMBRANDT® LSI Xp22.31-ISH detection assay C750K.XXXX.YY**

*All products are also available as stand-alone RUO products

**Available in 5 or 10 tests