

## ImmunoStrip® for Bt-Cry34Ab1



**Part Number:** STX 04500/0050

**Test Label:** NA  
**Test Format:** ImmunoStrip®  
**Capture Reagent:** Monoclonal  
**Detection Reagent:** Monoclonal

### Intended Use:

This ImmunoStrip® test is intended for seed quality purposes to determine the presence of the insect resistance trait (Bt-Cry34Ab1) in transgenic corn. The test is appropriate for testing individual seeds and leaves.

The Bt-Cry34Ab1 ImmunoStrip® has shown no cross-reaction with Bt-Cry1Ab, Bt-Cry1F, Bt-Cry3Bb1, Bt-Cry9C, EPSPS (Roundup Ready® events GA21 and NK603), or PAT/Liberty Link® transgenic proteins.

See product insert for additional details.

### Includes:

- ImmunoStrips®
- SEB4 sample buffer not included, see related products
- [User Guide](#)

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)