

ImmunoStrip® for Bt-Cry2Ab



Part Number: STX 05801/0050

Test Label: NA

Test Format: ImmunoComb® or ImmunoStrip®

Capture Reagent: Monoclonal

Detection Reagent: Monoclonal

Intended Use:

Especially designed for plant breeders, the Agdia Bt-Cry2A ImmunoStrip effectively detects the Bt-Cry2Ab protein from Monsanto's event MON89034

hybrid or inbred seed corn and leaves.

Once sample is prepared, test completion times are usually between 3 to 5 minutes.

The Bt-Cry2A ImmunoStrip® is applicable for use in the lab or the field, but special emphasis was placed on making a test useful for plant breeders under field conditions. All Agdia ImmunoStrips® are rugged, reliable and ready for field use.

Includes:

- ImmunoCombs®
- SEB4 sample buffer not included, see related products
- [User Guide](#)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com