

ImmunoStrip® for Bt-Cry2Ab & Bt-Cry1Ab/1Ac



Part Number: STX 06801/0050

Test Label: NA
Test Format: ImmunoStrip
Capture Reagent: NA
Detection Reagent: NA

Intended Use:

This ImmunoStrip® was developed for the cotton seed industry for those who need to test Bollgard® II varieties for purity. The ImmunoStrip® has two test lines; one for the Bt-Cry2A protein and one for the Bt-Cry1Ab/Ac protein. The ImmunoStrip® is suitable for testing single seeds and leaves. This test is also applicable to corn varieties being developed with the same proteins.

These ImmunoStrips are offered as single strips or in comb formats for 48 and 96 extraction plates.

Includes:

- ImmunoStrips
- [User Guide](#)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com