

ImmunoStrip® for Bt-Cry1F & Bt-Cry34Ab1



Part Number: STX 10900/0050
Test Label: NA
Test Format: ImmunoComb® or ImmunoStrip®
Capture Reagent: Monoclonal
Detection Reagent: Monoclonal

Intended Use:

This dual trait ImmunoStrip® test is intended for seed quality purposes to determine the presence or absence of the insect resistance traits (Bt-Cry1F and Bt-Cry34Ab1) HERCULEX® XTRA in transgenic corn.

The Bt-Cry1F and Bt-Cry34Ab1 dual trait ImmunoStrip® is appropriate for testing single seeds and leaves. Plant breeders in particular regularly use this dual-trait ImmunoStrip® for checking segregations.

In the field or in the lab, Agdia ImmunoStrips® give you clear results in about 5 minutes.

Includes:

- ImmunoComb® or ImmunoStrip®
- SEB4 sample buffer not included, see related products
- [User Guide](#)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com