

## ImmunoStrip® for Bt-Cry1F, Bt-Cry1Ac & CP4 EPSPS



**Part Number:** STX 13200/0050

### Intended Use:

This ImmunoStrip® test is intended for seed quality purposes to determine the presence of the Bt-Cry1F and Bt-Cry1Ac traits along with Roundup Ready® trait (CP4 EPSPS) in transgenic cotton seed and leaves.

This is a multi-analyte ImmunoStrip® that has three specific test lines. The Bt-Cry1F test line has shown no cross-reaction with other transgenic proteins in cotton seed and leaves including Bt-Cry1Ab, Bt-Cry1Ac, Bt-Cry2A, CP4 EPSPS (Roundup Ready®) or PAT/bar transgenic proteins.

The test line for Bt-Cry1Ac recognizes both Bt-Cry1Ac and Bt-Cry1Ab. It shows no cross reaction to Bt-Cry1F, Bt-Cry2A, CP4 EPSPS (Roundup Ready®) or PAT/bar transgenic proteins.

The CP4 EPSPS (Roundup Ready®) test line has shown no cross-reaction with Bt-Cry1Ab, Bt-Cry1Ac, Bt-Cry1F, Bt-Cry2A or PAT/bar transgenic proteins

### Includes:

- ImmunoStrips®
- SEB4 sample buffer not included, see related products
- [User Guide](#)

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)