


ImmunoStrip® for Hosta virus X (HVX)



Test Label: NA
Test Format: ImmunoStrip®
Capture Reagent: Polyclonal
Detection Reagent: Polyclonal

Intended Use:

The HVX ImmunoStrip® is an on-site tool to quickly identify plants infected with *Hosta virus X* (HVX). The HVX ImmunoStrip® can be used to test both leaves and roots of hosta plants.

The ImmunoStrip® test was validated using over 35 different varieties of Hosta. Several HVX isolates were tested including isolates from Dr. Lockhart and from several infected plants of commercial production sites.

The HVX ImmunoStrip® test is an ideal tool for on-site screening of plants. Use it in scouting regimes or to test plants before purchasing or accepting a shipment of plants.

Kit (ISK) Includes:

- ImmunoStrips®
- Sample bags containing SEB1 buffer
- [User Guide](#)

ImmunoStrips® Only (STX) Includes:

- ImmunoStrips®
- SEB1 extraction buffer purchased separately

Option

ImmunoStrip® for Hosta virus X (HVX), 5 strips and buffer filled bags
ISK 16600/0005

ImmunoStrip® for Hosta virus X (HVX), 25 strips and buffer filled bags
ISK 16600/0025

ImmunoStrips® only for Hosta virus X (HVX) ***BUFFER FILLED BAGS
NOT INCLUDED***
STX 16600/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com