

ImmunoStrip® for Tomato ringspot virus (ToRSV)



Test Label: NA
Test Format: ImmunoStrip®
Capture Reagent: Monoclonal
Detection Reagent: Monoclonal

Intended Use:

The ToRSV ImmunoStrip® is an on-site tool to quickly identify *Tomato ringspot virus* (ToRSV) in plant foliage.

ImmunoStrip® tests require no equipment or expertise to run. Results are obtained in as little as a few minutes making them perfect for use in the field or greenhouse.

Kit (ISK) Includes:

- ImmunoStrips®
- Sample bags containing SEB1 buffer
- [User Guide](#)

ImmunoStrips® Only (STX) Includes:

- ImmunoStrips®
- SEB1 extraction buffer purchased separately

Option

ImmunoStrip® for Tomato ringspot virus (ToRSV), 5 strips & buffer filled bags
ISK 22001/0005

ImmunoStrip® for Tomato ringspot virus (ToRSV), 25 strips & buffer filled bags
ISK 22001/0025

ImmunoStrips® only for Tomato ringspot virus (ToRSV) ***BUFFER FILLED BAGS NOT INCLUDED***
STX 22001/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com