

## ImmunoStrip® for Arabis mosaic virus (ArMV)

The *Arabis mosaic virus* (ArMV) ImmunoStrip® is used to detect the presence of ArMV, the causal agent of Raspberry yellow dwarf disease, in crops including grapevines, raspberry, and strawberry, as well as other fruit, vegetable, and ornamentals. ArMV is a member of the Nepovirus genus known for their transmission by nematodes. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

### Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

### Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

<input type="checkbox"/>	Option
<input type="checkbox"/>	ImmunoStrip® for Arabis mosaic virus (ArMV), 5 strips & buffer filled bags ISK 23203/0005
<input type="checkbox"/>	ImmunoStrip® for Arabis mosaic virus (ArMV), 25 strips & buffer filled bags ISK 23203/0025
<input type="checkbox"/>	ImmunoStrips® only for Arabis mosaic virus (ArMV) ***BUFFER FILLED BAGS NOT INCLUDED*** STX 23203/0025

### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)