

## ImmunoStrip® for AAD-1



**Part Number:** STX 23500/0050

**Test Label:** NA  
**Test Format:** ImmunoStrip®  
**Capture Reagent:** Monoclonal  
**Detection Reagent:** Monoclonal

### Intended Use:

This ImmunoStrip® test is intended for seed and leaf quality purposes to determine the presence or absence of the Enlist™ trait (AAD-1) that is resistant to 2,4-D herbicide in transgenic corn seed and leaf.

This ImmunoStrip® does not cross-react with other transgenic proteins in cotton including CP4 EPSPS, Bt-Cry1Ab/1Ac, Bt-Cry2A, Bt-Cry3A, mBt-Cry3A, Bt-Cry3Bb1, Bt-Cry34Ab1, Bt-Cry35Ab1, Bt-Cry1F, Bt-Cry9C, PAT/pat, GA21, GT21, CspB and Vip3A.

ImmunoStrip® tests require no expertise to run. Results are obtained in as little as a few minutes making them perfect for use in the field.

*Enlist™ and the Enlist logo are trademarks of The Dow Chemical Company. ("Dow") or an affiliated company of Dow.*

### Includes:

- ImmunoStrips®
- SEB4 sample extraction buffer not included, see related products
- [User Guide](#)

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)