

## ImmunoStrip® for Tomato spotted wilt virus (TSWV)



The *Tomato spotted wilt virus* (TSWV) ImmunoStrip® is used to detect the presence of TSWV in ornamental and vegetable crops. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

### Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

### Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

### Option

ImmunoStrip® for Tomato spotted wilt virus (TSWV), 5 strips & buffer filled bags  
ISK 39300/0005

ImmunoStrip® for Tomato spotted wilt virus (TSWV), 25 strips & buffer filled bags  
ISK 39300/0025

ImmunoStrips® only for Tomato spotted wilt virus (TSWV) \*\*\*BUFFER FILLED BAGS NOT INCLUDED\*\*\*  
STX 39300/0025

### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)