

ImmunoStrip® for Cucumber mosaic virus (CMV)

The *Cucumber mosaic virus* (CMV) ImmunoStrip® is used to detect the presence of CMV in ornamental and vegetable crops. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab. This test detects all known CMV isolates including all isolates belonging to subgroup I and subgroup II.

Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

| <input type="checkbox"/> | Option |
|--------------------------|--|
| <input type="checkbox"/> | ImmunoStrip® for Cucumber mosaic virus (CMV), 5 strips & buffer filled bags ISK 44501/0005 |
| <input type="checkbox"/> | ImmunoStrip® for Cucumber mosaic virus (CMV), 25 strips & buffer filled bags ISK 44501/0025 |
| <input type="checkbox"/> | ImmunoStrips® only for Cucumber mosaic virus (CMV) ***BUFFER FILLED BAGS NOT INCLUDED*** STX 44501/0025 |

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com