

ImmunoStrip® for eCry3.1Ab



Part Number: STX 49000/0050

Test Label: NA
Test Format: ImmunoStrip
Capture Reagent: Polyclonal
Detection Reagent: Monoclonal

Intended Use:

This ImmunoStrip® test is intended for seed quality purposes to determine the presence or absence of eCry3.1Ab protein in transgenic corn (Event 5307, Agrisure Duracade®).

This test system can be used to test individual corn seed and corn leaf or detect 1 transgenic eCry3.1Ab seed in 100 corn seeds (1.0 %).

The eCry3.1Ab ImmunoStrip® has shown no cross-reaction with other transgenic proteins in corn seed and leaf including GA21 (mEPSPS), Bt11 (Cry1Ab), mCry3A, or MIR162 (Vip3A).

Includes:

- ImmunoStrips®
- PBST sample buffer not included, see related products
- SEB4 extraction buffer powder not included, see related products
- [User Guide](#)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com