


ImmunoStrip® for Alternanthera mosaic virus (AltMV) / Papaya mosaic virus (PapMV)



Intended Use:

The AltMV / PapMV ImmunoStrip® is an on-site tool validated to detect Alternanthera mosaic virus (AltMV) in ornamental plant species.

This assay was tested against and accurately detects the following AltMV isolates: AltMV-SP, AltMV-PA, AltMV-LR, AltMV-GL, PV-1172, PV-0849, PV-0725. More than 80 healthy ornamental hosts were evaluated for diagnostic specificity and no host reactions (false positives) were observed.

This test can be used by growers, plant breeders, disease clinics, researchers, commercial laboratories, regulatory agencies and anyone testing for virus infection in developed plant, and plant part used for vegetative propagation.

Kit (ISK) Includes:

- ImmunoStrips®
- Sample bags containing SEB1 buffer
- [User Guide](#)

ImmunoStrips® Only (STX) Includes:

- ImmunoStrips®
- SEB1 extraction buffer purchased separately

Option

ImmunoStrip for AltMV/PapMV, 5 strips & buffer filled bags
ISK 53400/0005

ImmunoStrip for AltMV/PapMV, 25 strips & buffer filled bags
ISK 53400/0025

ImmunoStrip for AltMV/PapMV, 25 strips *** BUFFER FILLED BAGS NOT INCLUDED ***
STX 53400/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com